Dok.No / Doc. No	MDR.TL.15.15M(EN)
Yay.Trh / Pub. D.	22.04.2024
Rev.No / Rev. No	02
Rev.Trh. / Rev. D.	03.10.2024
Sayfa / Page	1 / 3

POWERBONE Instruction for Use of Low Viscosity (LV) Bone Cement**Classification**

- POWERBONE LV Bone Cement is a medical device (Class IIb), not a drug.

Material

- POWERBONE LV Bone Cement is a polymethylmethacrylate based radiopaque acrylic bone cement. Bone Cement is a reliable biomaterial due to its ease of use in clinical practice and its long survival rate, especially proven with cemented prostheses.

POWERBONE LV Product Components of Bone Cement:**Solid Phase**

- Poly (Methyl Methacrylate) (PMMA)
- Zirconium Dioxide (ZrO₂)
- Benzoyl Peroxide (BPO)

Liquid Phase

- Methyl Methacrylate(MMA)
- N,N-Dimethyl-p-toluidine (DMPT)
- Hydroquinone

General Description**Use and Features**

- POWERBONE LV Bone Cement is designed to be sent with a piece through vertebroplasty, kyphoplasty applications and cannulated pedicle screw. POWERBONE LV Bone Cement is not intended for use in paste form.
- Bone Cement is a polymer used in bone diseases that hardens and provides dispersed regional mechanical strength.
- The degradable polymer, which is formed as a result of the mixture of two components in liquid and solid phase, is collected into a special joining parts into the broken vertebra code and solidifies in the relevant area, providing mechanical strength.
- During injection, zirconium dioxide contained in POWERBONE LV Bone Cement provides injection performance in the appropriate area in the X-ray image.
- POWERBONE LV Bone Cement is manufactured and tested based on ISO 5833 International Standard at 23±1 °C.

Product General Features

- Biocompatible,
- Biomechanical support,
- Radiopaque,
- Easy implantation.

Indications

- Avascular necrosis
- Sickle cell anemia
- Fixation of pathological fractures where loss of bone substance or recalcitrance of the fracture renders more conventional procedures ineffective.

Contraindications

- PMMA Bone Cement is contraindicated in the presence of active or incompletely treated infection, at the site where the Bone Cement is to be applied. The use of PMMA is also contraindication for patients who:
 - Are pregnant or nursing
 - Are allergic to the any of the components of PMMA

- Have an active infectious arthritis of the joint or joints to be replaced or a history of such an infection
- Have a loss of musculature or have neuromuscular comprise in the affected limb this would render the procedure unjustifiable
- Have myasthenia gravis
- Have metabolic disorders which may impair bone formation
- Are hypotensive
- Have renal impairment
- Have congestive heart failure

Precautions

- Do not use this product after the expiration date printed on the package. This device may not be safe or effective beyond its expiration date. Follow the handling and mixing instructions to avoid contact dermatitis. Strict adherence to the instructions for mixing the powder and liquid components may reduce the incidence of this complication. Adequately ventilate the operating room to eliminate as much monomer vapor as possible. The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted Bone Cement has been reported. Since soft contact lenses are quite permeable, they should not be worn in an operating room where methyl methacrylate is being mixed. Bone Cement is only intended for use by specialized orthopedic doctors, neurologist and neurosurgeons.

Possible Adverse Effects

- Serious adverse events, some with fatal outcome, associated with the use of acrylic Bone Cement include myocardial infarction, cardiac arrest, acute stroke and pulmonary embolism.
- The most frequent adverse reactions reported with acrylic Bone Cement are a transitory fall in blood pressure, thrombophlebitis, hemorrhage and hematoma, loosening or displacement of the prosthesis, superficial or deep wound infection, trochanteric bursitis, and short-term cardiac conduction irregularities. Other reported adverse reactions include heterotopic new bone formation and trochanteric separation.
- Other reported adverse events for acrylic Bone Cement include pyrexia due to an allergy to the Bone Cement, hematuria, dysuria, bladder fistula, delayed sciatic nerve entrapment due to extrusion of the Bone Cement beyond the region of its intended application, and adhesions and stricture of the ileum due to the heat released during polymerization.

Warnings

- Monitor patients carefully for any change in blood pressure during and immediately following the application of Bone Cement. Adverse patient reactions affecting the cardiovascular system have been associated with the use of Bone Cement s. Hypotensive reactions have occurred between 10 and 165 seconds following application of Bone Cement; they have lasted from 30 seconds to 5 or more minutes. Some have progressed to cardiac arrest. Patients should be monitored carefully for any change in blood pressure during and immediately following the application of Bone Cement. Follow the handling, mixing and canal preparation instructions carefully.
- Caution should be exercised during the mixing of the two components to prevent excessive exposure to the concentrated monomer vapors, which may produce irritation of the respiratory tract, eyes, and possibly the liver. Personnel wearing contact lenses should not be near or involved in mixing this Bone Cement.

Hazırlayan / Prepared By	Kontrol eden / Controlled by	Onaylayan/ Approved By
KALİTE GÜVENCE PERSONELİ QUALITY ASSURANCE PERSONNEL	KALİTE GÜVENCE SORUMLUSU QUALITY ASSURANCE RESPONSIBLE	KALİTE YÖNETİM TEMİSİL CİSİ QUALITY MANAGEMENT REPRESENTATIVE

Dok.No / Doc. No	MDR.TL.15.15M(EN)
Yay.Trh / Pub. D.	22.04.2024
Rev.No / Rev. No	02
Rev.Trh. / Rev. D.	03.10.2024
Sayfa / Page	2 / 3


- The hardening of Bone Cement is an exothermic reaction and it releases heat. The heat generated as a result of this reaction can damage the applied bone tissue.
- When working with the monomer or the cement, gloves must be worn to ensure adequate protection against the penetration of the monomer (methyl methacrylate) into the skin. PVP (three-layer polyethylene, ethylene-vinyl alcohol copolymer, polyethylene) and Viton/Butyl gloves have proved to provide good protection over an extended period. Putting on two pairs of gloves – one pair of polyethylene surgical gloves over a pair of standard latex surgical gloves has also proved to offer adequate protection. The use of latex or polystyrene-butadiene gloves, however, must be avoided. Please request the confirmation of your glove supplier whether the respective gloves are suitable for the use of this cement.
- Avoid over pressurization of the Bone Cement because this may lead to extrusion of the Bone Cement beyond the site of its intended application and damage to the surrounding tissues.
- The safety of the Bone Cement in pregnant women or in children has not been established. Bone Cement may adversely affect bone growth and fetal health.

Dosage

- The amount of Bone Cement required depends on the patient's anatomy and on the implant used. One or more complete units (the contents of one pouch and one ampoule) must always be mixed together. In practice, more than four dose units per implantation are only rarely used.

Application Notes on Use

- Before using POWERBONE LV Bone Cement the surgeon must be thoroughly familiar with its properties, handling and use, and must have read the user manual. Because the handling and curing characteristics of this Bone Cement vary with temperature, humidity, and mixing technique.
- POWERBONE LV Bone cement mixing, syringe filling, working time and hardening times are given in the chart below.

- Before using POWERBONE LV Bone Cement for the first time, a test mixing should be carried out to become familiar with the characteristics of BONEGRAFT Bone Cement.
- POWERBONE LV Bone Cement and equipments to be used should be stored in the application environment at least 2 hours before use.
- After the product is removed from its protective packaging, the outer and inner sterilization roll package of the powder phase and the sterile outer package of the liquid phase (Pet-g, Tyvek) should be opened by the circulating nurse, maintaining sterility, and the ampoule package should be opened for mixing. The ampoule is opened by breaking its neck. Do not open the ampoule over the mixing device to prevent contamination of the cement with glass fragments.

Opening under sterile conditions:

- The opening flaps from the bag top assist to detach the foil from the paper. In order to grip as much of the opening flaps as possible, the side of the paper / foil should be kept between thumb, index and middle finger. Please use the whole thumb surface to grip foil and paper side and take off each side evenly.

Mixing Time

First pour all the liquid from the ampoule into a sterile mixing device and add all the powder. Stir carefully with a sterile mixing rod until a homogeneous mass is obtained. The mixture should be mixed for 30 sec, irrespective of the room temperature. The temperature/time diagrams must be consulted. Once the mixing process completed initial dough is reached and syringe filling for vertebroplasty/kyphoplasty applications can be conducted.

Syringe Filling

- During this stage syringe filling for vertebroplasty/kyphoplasty applications can be conducted.

Working Procedure

- After careful preparation for Kyphoplasty or Vertebroplasty applications, POWERBONE LV Bone Cement can be applied using a cement syringe or an alternative injection delivery system. Method of application is given in the instructions for the system used and must be consulted. The working time and the rate of polymerization are strongly dependent on ambient temperature and on temperature of the components. The hardening time is shorter at high temperatures and longer at low temperatures. Viscosity increases as polymerization progresses, i.e. as the cement hardens.

Storage and Preservation

- Do not store above 25°C (77°F). Store in a DRY PLACE AT ROOM TEMPERATURE. Optimal storage conditions: +5°C - +25°C and normal relative humidity (30-65%).

Shelf Life/ Sterility

- The expiration date is printed on the outer carton, the protective aluminum packaging and the inner pouch. The expiration date is 3 years and printed on the label. **DO NOT USE AFTER THE EXPIRATION DATE.**
- POWERBONE LV Bone Cement must not be used after the expiration date. The contents of unused but opened or damaged packs must not be resterilized and are therefore to be discarded.
- POWERBONE LV Bone Cement is sterilized with ethylene oxide gas and must not be resterilized.

Disposal


- Single components of the Bone Cement, cured solid material as well as (uncleaned) packaging material must be disposed by following the regulations of the local authorities. Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

Further Information:

- For further information on the product and its uses, please contact BONEGRAFT Biyolojik Malzemeler San. Tic. ve A.Ş.. The address is printed on this information sheet.

Pictogram For Product

Pictogram		Explanation
Hazırlayan / Prepared By KALİTE GÜVENCE PERSONELİ QUALITY ASSURANCE PERSONNEL	Kontrol eden / Controlled by KALİTE GÜVENCE SORUMLUSU QUALITY ASSURANCE RESPONSIBLE	Onaylayan/ Approved By KALİTE YÖNETİM TEMSİLCİSİ QUALITY MANAGEMENT REPRESENTATIVE

	Expire Date
	Sign of Conformity
	Do Not Use If Damaged
	Sterilized by EO for powder
	Sterilized by Aseptic for Liquid
	Consult Instruction for use
	Temperature Limit
	Humidity Limit
	Manufacturer
	Lot number / Batch Code
	Do Not Reuse
	Do not Re-Sterilize
	Catalog Number
	Date of Manufacture Manufacturer Country
	Warnings
	Keep Away From Sunlight
	Keep Dry
	Inflammable Matter

	Medical Device
	Double Sterile Barrier System
	Unique Identifier Information
	QR

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Hazırlayan / Prepared By KALİTE GÜVENÇE PERSONELİ QUALITY ASSURANCE PERSONNEL	Kontrol eden / Controlled by KALİTE GÜVENÇE SORUMLUSU QUALITY ASSURANCE RESPONSIBLE	Onaylayan/ Approved By KALİTE YÖNETİM TEMSİLCİSİ QUALITY MANAGEMENT REPRESENTATIVE
--	--	---

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