Data 16-06-2022

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
Signature

 Mr. Vivek Kumar Shahi
 Regulatory Affairs Assistant
 Aegis Lifesciences Pvt. Ltd.
 Place: Ahmedabad, India

**Technical Data sheet for Surgi-ORC®**

Company details	Name of Manufacturer Address Telephone number Fax E-mail Address	Aegis Lifesciences Private Limited 215 - 216, Mahagujarat Industrial Estate, At: Moraiya, PO: Changodar, TAL: Sanand, Ahmedabad – 382213 Gujarat, India. +91 9726903131 +91 9727123131 admin@aegis-lifesciences.com
Products description and indication	Surgi-ORC® absorbable haemostat is a sterile absorbable knitted fabric prepared by controlled oxidation of regenerated cellulose. The fabric is white to pale yellow in color and has a faint, caramel-like aroma. It is stable at room temperature. A slight discoloration may occur with age, but this does not affect its performance. The Surgi-ORC® haemostats are available as Original/Standard, Knit, Fibril and Non-Woven/SNOW, manufactured with the same Oxidized Regenerated Cellulose. Surgi-ORC® meets all the requirements of Oxidized Regenerated Cellulose established by United States Pharmacopoeia (USP). The Surgi-ORC® haemostat is used adjunctively in various surgical procedures to assist when control of bleeding from capillary, venous and small arteriolar vessels, by pressure, ligature and other conventional procedures are either ineffective or impractical. The Surgi-ORC® haemostat can be cut to size in endoscopic procedures.	
Directive Class (under Directive 93/42/EEC)	III	Rule of classification: Rule 8 (Annex IX, MDD 93/42/EEC)
Mode of action	When Surgi-ORC® haemostat contacts with blood, it forms a brownish or black gelatinous mass, which aids in formation of clot, thereby serving as a haemostatic adjunct in the control of local haemorrhage. When Surgi-ORC® haemostat is used properly in minimal amounts, it is absorbed from the sites of implantation without tissue reaction. Absorption depends upon several factors including the amount used, degree of saturation with blood and tissue bed. In addition to its local haemostatic properties, Surgi-ORC® haemostat is bactericidal in vitro against a wide range of gram positive and gram negative organisms including aerobes and anaerobes. However, Surgi-ORC® haemostat is not an alternative to the systemically applied therapeutic or prophylactic antimicrobial agents.	
Technical Specifications	Product is sterile, white to pale yellow in color and has a faint, caramel-like aroma, biocompatible, acidic in nature with a pH value of 1-3 and Haemostasis Time of 3 to 4 minutes.	
Size/type:	Please refer to table below.	

Sr.No.	Model Variant	Product ID	Dimensions	Number of pieces per packaging
1.	Original /Standard	SOO-502	0.5 in x 2 in (1.25 cm x 5.0 cm)	12 Pcs / Box
2.		SOO-0203	2 in x 3 in (5.1 cm x 7.6 cm)	12 Pcs / Box
3.		SOO-0214	2 in x 14 in (5.0 cm x 35 cm)	12 Pcs / Box
4.		SOO-0408	4 in x 8 in (10 cm x 20 cm)	12 Pcs / Box
5.	Fibril	SOF-0102	1 in X 2 in (2.5 cm X 5.1 cm)	12 Pcs / Box
6.		SOF-0204	2 in X 4 in (5.1 cm X 10.2 cm)	12 Pcs / Box
7.		SOF-0404	4 in X 4 in (10.2 cm X 10.2 cm)	12 Pcs / Box
8.	Non- Woven/SNOW	SON-0102	1 in x 2 in (2.5cm x 5.1 cm)	12 Pcs / Box
9.		SON-0204	2 in x 4 in (5.1 cm x 10.2 cm)	12 Pcs / Box
10.		SON-0404	4 in x 4 in (10.2 cm x10.2 cm)	12 Pcs / Box
Composition / ingredients		Oxidized Regenerated Cellulose		
Clinical safety		Product has history of safe use. The US Pharmacopoeia contains monograph for Absorbable Gelatin Sponge. Conformity with the requirements concerning the characteristics and performances referred to clinical evidence guideline with European Commissions Guideline MEDDEV 2.7/1, Rev.4.under the normal conditions of use of the device and the evaluation of the undesirable side-effects has been achieved in conformity with the requirement of MDD 93/42/EEC.		
Biocompatibility		The device is biocompatible. Biocompatibility tests is carried out as per international standard ISO 10993 Series.		
Sterilization		The product, Surgi-ORC® is supplied in sterile form. Sterilization is carried out by gamma irradiation. It is determined that Gamma irradiation have no effect on the materials used in the product, its components and packaging.		
Packaging description		Surgi-ORC® is packed in a complete sterile barrier packaging system having high quality/medical grade quality packing material for all types of packing stages.		
Label		As per ISO 15223-1:2016		
Instructions for use		See Attached IFU		
Precautions for handling or storage		Store the product in its original packaging in a clean and dry room at the temperature not more than 30°C. Close the outer package of the product immediately after use. Do not refrigerate or freeze.		
Traceability		GS1 Barcode		

Mr. Bhavin Trivedi
Regulatory Affairs Executive
Aegis Lifesciences Pvt. Ltd.
Ahmedabad, India.
Date: 22-10-2021